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The screenshot displays the Nikon NIS-Elements AR software interface. The main window shows the 'ND Acquisition Wizard' dialog box, which is used to configure time-lapse acquisition parameters. The wizard includes sections for 'Acquisition Parameters', 'Timelapse', 'Multipoint', and 'Z-Stack'. The 'Timelapse' section is currently selected, showing a table with columns for 'Name', 'Phase', and 'Interval'. The table contains one entry: '#1' with 'Fixed Duration' phase and '10 sec' interval. The 'Run' button is highlighted, indicating the start of the acquisition process.

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The screenshot displays the Nikon NIS-Elements AR software interface, specifically the 'Scanning Wizard' dialog box for large image stitching. The interface is divided into several panels and sections:

- Top Panel:** Contains icons for Live, Capture, Find, ND Acquire, and a 'Large Image' icon (highlighted by callout 11.1).
- Left Panel:** Shows ROI Statistics and Intensity Profile graphs.
- Scanning Wizard Dialog:**
 - 11.2:** 'Use Z for Focus Surface' checkbox is checked.
 - 11.3:** 'Save file(s) into folder' is selected, with the path 'C:\Program Files\NIS-Elements\Images'.
 - 11.4:** 'Shape' is set to 'M x N'.
 - 11.5:** 'Fields' are set to 4 x 2, and 'Overlap' is 10%.
 - 11.6:** 'Stitching via' is set to 'Blending'.
 - 11.7:** 'Run' button is highlighted.
- Right Panel:** Shows various acquisition parameters like Speed, Quality, Averaging, Dwell Time, Channel mode, and Filters.
- Bottom Panel:** Displays system information and a status bar with coordinates (X, Y, Z).

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