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Compact Advanced Wide

File Edit Acquire Calibration Image ROI Binary Measure Reference Macro View Devices Window Applications Assays HCA/JOBS Addons NIS.ai Deconvolution Help

7 564

G: 1.00

0 10000 20000 30000 40000 50000 60000

76 (6, 1252)

G: 1.00

0 10000 20000 30000 40000 50000 60000

20 (0, 920)

Sample Navigation

Stage Document

Full ROI Scan Overview Select Holder

T12-S HU Universal Holder: Slide Focus Position Set Reset

No data Scan for Sample Overview

Focus 1800 1600 1200 500

Escape

ROI Statistics Intensity Profile Measurements

# Length [px] EqDiameter [px] Count

Mean - - -

StDev - - -

Min - - -

Max - - -

Count 0 0 0

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Capture and Store Pad

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20231205\_multi-chns.nd2\* X Captured 1 : Denoised... X Captured 6.nd2\*

AX

11 91%

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AX GREEN RED 0.04 um/px 2x16bit: 2048 x 2048 pixels [N/A]

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Live Capture Find AutoSignal Autofocus Short Autofocus Long ND Acquire Large Image

DIA D-LEDE AX AX NSPARC 光路切换 Lightpath

Experiment: Search Sample Settings

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The screenshot displays the Nikon NIS-Elements AR software interface. The main window shows a live view of a sample with multiple channels (DAPI, EGFP, RFP). Overlaid on this are several windows and callouts:

- ND Acquisition Wizard:** A central window with tabs for Acquisition Parameters and Timelapse. The Timelapse tab is active, showing options for Timelapse, Multipoint, Z-stack, and Multichannel. The 'Filename' field is set to 'Experiment'.
- Callout 1:** Points to the 'PFS' (Perfect Focus System) button in the bottom right corner of the main window.
- Callout 2:** Points to the 'ND Acquisition' button in the top right corner of the main window.
- Callout 3:** Points to the 'Timelapse' checkbox in the ND Acquisition Wizard.
- Callout 4:** Points to the 'Timelapse' tab in the ND Acquisition Wizard.
- Callout 5:** Points to the 'Save file(s) into folder' radio button in the ND Acquisition Wizard.
- Callout 6:** Points to the 'Time Sequence Name' field and the table below it in the ND Acquisition Wizard.
- Callout 7:** Points to the 'Run' button in the bottom right corner of the ND Acquisition Wizard.

**Time Sequence Name Table:**

Name	Phase	Interval
#1	Fixed Duration	3 min
		Run Every
		1 sec

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The screenshot displays the Nikon NIS-Elements AR software interface. The main window shows a live view of a sample with green and red fluorescence. Overlaid on this are several windows and panels:

- ND Acquisition Wizard:** A central dialog box with tabs for Acquisition Parameters, Timelapse, and Multipoint. The Timelapse and Multipoint options are checked. The Multipoint section shows a list of predefined points with their X, Y, and Z coordinates.
- Acquisition Parameters:** A panel on the left showing various acquisition settings like exposure time, gain, and filter.
- Timelapse:** A panel on the left showing settings for time-lapse imaging.
- Multipoint:** A panel on the left showing settings for multi-point imaging.
- ROI Statistics:** A panel at the bottom left showing statistics for the region of interest.
- Measurements:** A panel at the bottom right showing measurement results.

Yellow callout boxes provide numbered instructions for setting up the acquisition:

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